

BVT02 Bladder Scanner

- Features
- Multi-language customized
- ② 3.5' color LCD screen, good interface style
- 3 Three patients figure mode (Standard / Obesity/ Children)
- ④ Scan time <4S</p>
- ⑤ Machine beeps when the measurement exceeds the capacity of the alarm set value
- Select different image contrast, can improve the measurement accuracy
- Probe preset time automatically scan within specified time
- Work time≥9 hours

Application & Effects

Home health

Control the micturition time of patients by measuring post-void residual volume, avoiding the irrevocable damage of bladder function caused by long-term open catheterization

Geriatrics

Test urine volume can evaluate urinary retention , uracratia and other urinary system diseases

Rehabilitation care

Prevent urinary tract infection caused by unnecessary catheterization and reduces patient pain

Urology

Evaluate the therapeutic effect of certain drugs for urinary system disease by continuously measuring post-void residual volume




KALXIN (





Standard configuration

Main unit 1pc 2.5MHz 3D probe 1pc Power adapter 1pc Internal batter 1pc



calibrmtion tool



ADD:Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.

Zip Code: 221004

Tel: +86-516-87732932 / 87733758

Fax: +86-516-87732932 / 87792848

http://www.kaixinbvt.com Http://www.kxele.com E-mail:info@kxele.com / service@kxele.com

